# Determination of critical micellar concentration by NMR spectroscopy

Laboratory Practice,

700 MHz NMR spectrometer, -1.108 lab

Practical: Dr. Bodor Andrea

abodor@caesar.elte.hu



## Goal

Critical micellar concentration of an anionic surfactant, sodium dodecyl sulphate (SDS) is going to be determined using diffusion NMR measurements.



## Micelle formation

![](_page_2_Figure_1.jpeg)

Critical micellar concentration (CMC)

The minimum required concentration of the surfactant to initiate micelle formation. Above this concentration the amphiphilic molecules form micelles.

![](_page_2_Figure_4.jpeg)

Hagn, F.; Etzkorn, M.; Raschle, T.; Wagner, G. J Am Chem Soc 2013, 135, 1919–1925.

### Determination of cmc

![](_page_3_Figure_1.jpeg)

Today: cmc is going to be determined from diffusion coefficients. Diffusion coefficient determination using NMR.

### **Stokes-Einstein equation**

![](_page_4_Figure_1.jpeg)

- D-diffusion coefficient
- k Boltzmann constant (1.38·10<sup>-23</sup> J/K)
- T-temperature
- $\eta$  viscosity
- $r_{\rm H}$  hydrodynamic radius

The diffusion coefficient is inversly proportional to the size of the particle.

![](_page_4_Figure_8.jpeg)

The measured diffusion coefficient is the mole fraction weighted average of free and in micelle surfactants.

$$D_{\text{measured}} = x_{\text{free}} D_{\text{free}} + x_{\text{mic}} D_{\text{mic}}$$

How does the diffusion coefficient change, when the total surfactant concentration increases?

## **PGSE-Pulsed Gradient Spin Echo**

#### Schematic representation of the pulse program

![](_page_5_Figure_2.jpeg)

We used **stebpgp1s19** pulse program for the DOSY measurements DOSY: **D**iffusion **O**rdered **S**pectroscop**Y** 

### **DOSY** measurements

![](_page_6_Figure_1.jpeg)

without diffusion

Determination of Diffusion coefficients: Stejskal-Tanner equation

$$I = I_0 \exp[-D\gamma^2 \delta^2 G^2 \left(\Delta - \frac{\delta}{3}\right)]$$

![](_page_7_Figure_2.jpeg)

Parameters to set up:  $\delta$ ,  $\Delta$ In one experiment we change: Gradient strength From a non-linear fitting we get: **Diffusion coefficient**.

## The steps of our present practical

- 1. Sample preparation
- 2. NMR measurements
  - 1. Lock
  - 2. Atma proton
  - 3. Shim
  - 4. Acquisition of 1D <sup>1</sup>H spectrum
  - 5. DOSY measurement
- 3. Evaluation:
  - 1. Fitting the Stejskal-Tanner equation (1D evaluation)
  - 2. EDDOSY (2D evaluation)

## Sample preparation

Reagents:

- 0.1 M SDS stock solution
- Distilled water
- D<sub>2</sub>O

![](_page_9_Figure_5.jpeg)

Prepare 5 diluted SDS solutions (Find the concentrations given in the table below)! Each sample has a total volume of 600  $\mu$ l and each contains 10 % D<sub>2</sub>O!

Calculate the amount of SDS, H<sub>2</sub>O and D<sub>2</sub>O in each sample!

c <sub>sDS</sub> / mM	V <sub>sDS</sub> / μl	V <sub>H2O</sub> / μΙ	V <sub>D2O</sub> / μΙ
3			
6			
9			
12			
15			

### NMR measurements

Software used for data processing and evaluation:

TopSpin 3.6.2

Download (free for students): https://www.bruker.com/service/support-upgrades/softwaredownloads/nmr/free-topspin-processing/nmr-topspin-license-foracademia.html

> Download **version 3.6** not version 4!!! Choose "Data processing only" ("Setup type")

#### TopSpin 3.6.2

#### A brief tutorial:

#### https://www.youtube.com/watch?v=FocoABJ2rvw

![](_page_11_Figure_3.jpeg)

Write commands here

Adding a folder: Right click on Spectrum files  $\rightarrow$  Add new data dir. Load a Spectrum: drag-and-drop to the spectrum window In the virtual laboratory practice you will be guided through the steps of recording <sup>1</sup>H 1D and diffusion NMR spectra (slides 14-25).

You will get a detailed hands-on on how to evaluate the spectra (slides 27 -40) and you have to perform this data analysis based on the spectra you receive in a separate file.

![](_page_12_Picture_2.jpeg)

### Measurement I – new dataset

edc: creating a new dataset with copying the parameters of a previous one

🖕 Create New Dataset - new	×		Name of folder:
Prepare for a new experiment by creating a new data initializing its NMR parameters according to the select For multi-receiver experiments several datasets are c Please define the number of receivers in the Options	set and ted experiment type. ;reated.		<sup>/</sup> Practice_Pr_[YYMMDD]
NAME	Gyakorlat_Pr_190923	<b>`</b> /	<ul> <li>Number of experiment in folder</li> </ul>
EXPNO	12		
PROCNO	1		Following datasets should be created:
Use current parameters			• <b>zgpr</b> : 1D proton with water
O Experiment	Select		
<ul> <li>Options</li> </ul>			suppression, only for setting O1
Set solvent	H2O+D2O ~		and pr (only 1)
C Execute 'getprosol'			<ul> <li>zgesgp: 1D proton with better</li> </ul>
◯ Keep parameters	P 1, O1, PLW 1 V Change		water suppression (1 for each
DIR	C:\Egyetem\ELTE\PhD\Oktatas\Nagymuszeres_DOSY ~		sample)
Show new dataset in new window			• <b>stebpgp1s19</b> : DOSY experiment
Number of additional datasets: (1,2,16)	1		with water suppression (1 for
		-1	each sample)
			Title: add sample properties here
	OK Cancel More Info Help		

## Measurement II – preparations

- 1. lock: monitoring the solvent's deuterium signal in order to compensate the temporal variations in the magnetic field  $\rightarrow$  we have to set the solvent
- 2. atma proton: Tuning the RF coil to the <sup>1</sup>H resonance frequency
- shim: ensures the spatial homogenity of the magnetic field → we compensate the inhomogenities of the 14.6 T magnetic field with shim-coils (procedure: Z, Z<sup>2</sup>, Z, Z<sup>2</sup>... then Z<sup>3</sup>, X, Y)

![](_page_14_Figure_4.jpeg)

### Measurement III – p1, O1, 1D <sup>1</sup>H spectrum

- *zg*: starting the measurement ← start the previously created **zgpr**
- *pulsecal*: calculates the length of the 90 ° hard pulse on the proton channel

![](_page_15_Figure_3.jpeg)

• *O1*: Transmitter frequency offset in Hz. We should set it to the maximum of the water signal to achieve maximum efficiency of water suppression.

![](_page_15_Figure_5.jpeg)

Good O1 - effective water suppression: residual water signal is small and distorted

### Measurement IV – 1D proton

- Create a **zgesgp** dataset
- getprosol 1H [p1] -12.55:
  - Insert p1 value determined with *pulsecal* to [p1]
  - a command that <u>sets *p1*</u> and the power of shaped pulses
- *O1* Set the O1 value determined from the **zgpr**
- zg start experiment  $\rightarrow$  what do we get?
- *efp* Exponential multiplication + <u>Fourier</u>-<u>transformation</u> + phase correction → what do we get after this command?

![](_page_16_Figure_8.jpeg)

#### Importance of water suppression

![](_page_17_Figure_1.jpeg)

#### Assignment of 1D <sup>1</sup>H spectra 0\_0 5 11 3 9 12 2 10 8 6 [Iel] 298K, SDS 3mM 4. 9 9 5. 3. 1. 40

water ~4,8 ppm 5 4 3 2 1 0 0 [ppm] Signals belonging to SDS We detect 5 signal groups To which SDS protons do the signals belong?

### Assignment should be included in the report

Number of signal (see previous slide)	Chemical shift (give accurate values in the report)	Assignment	Multiplicity	Relative integral
	~4 <i>,</i> 8 ppm	H <sub>2</sub> O	-	-
1.	~3 <i>,</i> 9 ppm			
2.	~1,6 ppm			
3.	~1,3 ppm			
4.	~1,2 ppm			
5.	~0,8 ppm			

### Diffusion measurements

Following parameters should be set:

General: p1, O1, SW etc.

<u>DOSY parameters</u>:  $\delta$ ,  $\Delta$ , gradient strength (G) is changed according to what function and in how many steps?

Setting  $\delta$ :

command: P30: we must give the half of  $\delta$  in  $\mu$ s!

Setting  $\Delta$ :

command: D20: we must give the  $\Delta$  value in s!

Changing gradient strength: See "Starting DOSY experiment" slide

![](_page_20_Picture_9.jpeg)

Pay attention to the units! When setting the parameters, half of  $\delta$  is in  $\mu$ s,  $\Delta$  is in s! When evaluating DOSY data,  $\delta$  and  $\Delta$  values must be given in ms!

### **DOSY** parameters

• Recommended values for  $\delta$  and  $\Delta$ :

- Small molecules:  $\delta$  = 2 ms,  $\Delta$  = 75, 100 ms P30 = 1000, D20 = 0.75, 0.1

- Larger molecules (10 kDa <):  $\delta$  = 4 ms,  $\Delta$  = 200 ms P30 = 2000, D20 = 0.2

For SDS:  $\delta = 2 \text{ ms}$ ,  $\Delta = 150 \text{ ms}$ 

## Starting DOSY experiment

#### (not needed at home)

🕹 Bruker TopSpin 3.6.2 on FANNI-PC as FANNI [Academic License]								_ 0 ×
<u>Start</u> <u>P</u> rocess A <u>n</u> alyse P <u>u</u> blish	<u>V</u> iew	<u>M</u> anage	2					1
Create Datase	et 📃 Find (	Dataset	Open Dat	aset 「 Paste	Dataset Read	Pars.		
		Hz	[H_x]					
		ppḿ L→ IIII						
Browser Last50 Groups	<b>1</b> Gyakorlat_Pr	190923 2 1	C:\Users\FANNI\Doc	uments\NMR\spektrun	ı			
C:\Bruker\TopSpin3.5pl7\examdata	Spectrum F	ProcPars A	cquPars Title P	ulseProg Peaks I	ntegrals Sample Struc	ture Plot Fid		
Spektrumok	298K, SDS	3mM						[i=
🐵 🚹 gyakorlat								
								-2
								-6
								E
								- 00
								- 0
					1			=
							··· <b>/</b> ····	- 4
								-
						/	_ال	
Structure			÷	<u>-</u>				
No structure available.	10		8	6	4	2		U [ppm]
		-						

Type *dosy* command here

A new panel will appear, where we can set the parameters of how to change the gradient strength.

🖕 dosy	×
8	Define the gradient ramp for shape SMSQ10.100 (Integral = 0.900)
	Close

Press ENTER here.

#### We set the following parameters:

1. First gradient amplitude, percentage of the maximum gradient strength

🤹 dosy 📃 🔀	🧅 dosy
Enter first gradient amplitude:	Enter final gradient a
2	98
OK Cancel	

3. In how many points do we change gradient strength

![](_page_23_Picture_4.jpeg)

2. Final gradient amplitude

🖕 dosy	×
Enter final gradient amplitude:	
98	
	OK Cancel

4. According to what function do we modify gradient strength

🧅 dosy	×
ramp type (I q e {linear/squared/	exponential} ):
	OK Cancel

During practice we modify gradient strength linearly in 8 steps, between 5 and 95%.

### A FID of a DOSY experiment

![](_page_24_Figure_1.jpeg)

8 points, 8 FIDs are seen with decreasing intensities

#### Evaluation of the diffusion measurements

Pseudo-2D dataset: Fourier-transform with the *xf2* command

![](_page_25_Picture_2.jpeg)

2D view

<sup>3</sup>D view

#### Evaluation based on fitting the Tanner equation

![](_page_26_Figure_1.jpeg)

T1/T2 row : click on <u>Fid</u>

icon: what appears

🖕 E	xtract a row from 2d data	×
?	Fid or Spectrum must be extracted From the 2d relaxation data. This row should correspond to an experiment with the maximum or minimum delay time All further data preparation will be done in respect to this row.	k
	FID Spectrum Cancel	
4		_
Spe	ectrum slice must be extracted From the 2d relaxation data.	,

All further data preparation will be done in respect to this spectrum.

Slice Number = 1	]	
	ОК	Cancel

Choose "Spectrum": this means you will pick one 1D spectrum from the pseudo 2D data

Choose the first spectrum (recorded with the lowest gradient intensity): 1 Then press OK.

#### The corresponding 1D <sup>1</sup>H spectrum appears. Zoom on the aliphatic region of the SDS peaks.

![](_page_28_Picture_1.jpeg)

Click now on the <u>A Peaks/Ranges</u> menu, and from the appearing window choose the "Manual integration" (next button will be ["Prepare relaxation data"] OK).

choose the integration (second icon from left), if this is clicked, then you can integrate the selected regions (click and drag on the spectrum) and in red the integral values will

![](_page_29_Figure_2.jpeg)

Click on the Relaxation icon. Don't worry an error-message looking-like text (complaining about a vd file) will pop-up, just click on it, and ignore. And now you will see the following table.

The first 9 entries correspond to the already acquired spectra, and deal with how to integrate them. We want to integrate all 8 of them, starting from the first slice and moving in one increments.

The following 5 entries will help the program define what type of equation to fit on the integrated intensities. As we varied the gradient strength based on the predefined difflist these should be set (see red boxes). Now we chose the Stejksal Tanner equation, and we have to tell which were the experimentally used diffusion parameters in the coming 4 entries (out of which only 3 are active – shown also by red boxes): these are the GAMMA, Little delta and Big delta.

![](_page_30_Figure_3.jpeg)

![](_page_31_Figure_0.jpeg)

### Evaluation 1 of the diffusion parameters

Data analysis presented on slides 27-33 has to be performed for each solution, pick the corresponding measurement (in TopSpin 1D spectrum measurement files are shown in black the ones referring to diffusion data are pseudo 2D and are shown in blue).

In the end you will have a *D* value for each SDS concentration, and this is the graph you have to make. The intersection of the two fitted linear equations will give you the cmc.

![](_page_32_Figure_3.jpeg)

Note how given peaks vary in the 1D <sup>1</sup>H spectra! (chemical shift variation). Which environments are mostly affected? Can you give an explanation? (You can overlay the 1D proton spectra with the 3 icon)

### Evaluation 2 of the diffusion parameters eddosy command

In the end you will get an NMR "cromatogram" looking like 2D table, where the vertical axis is the log of the diffusion coefficient (log*D*) and the horizontal axis is the <sup>1</sup>H 1D spectrum. In case of multicomponent system each molecule will show up as a horizontal line – its <sup>1</sup>H chemical shifts - with the corresponding log*D*. We show an example how can you determine the different components from a mixture:

![](_page_33_Figure_2.jpeg)

http://chem.ch.huji.ac.il/nmr/techniques/other/diff/diff.html

## eddosy evaluation

- 1. Recall the original diffusion measurement (close all other windows opened in TopSpin)
- 2. xf2 command (this is the processing, performs Fourier transformation for each spectra)
- 3. eddosy command

The following window appears, click on OK button

![](_page_34_Picture_5.jpeg)

#### The appearing window:

1 Gyakorlat_Pr_190923 3 1 C:\Users\FANNI\Documents\NMR\spektrum					
Spectrum	ProcPars AcquPars	Title PulseProg Peaks	Integrals Sample Structure Plot Fid		
PG	L I; ▶ U V (	<i>e</i>			
General First	General				
Second	Method	exponential -	Processing method	=	
Third	ExpVar	Gradient -	Variable parameter		
Baseline	Xlist	difflist	Variable parameter values file name		
Contin	Nstart	0	Start of input points		
	Ndata	8	Number of input points (TD)		
	Maxiter	100	Maximum number of iterations		
	EPS	1	Tolerance		
	Nexp	1	Number of components to fit		
	Noise	24351.00	Noise level (S_DEV)		
	PC	4	Noise sensitivity factor		
	SpiSup	1	Spike suppression factor		
	F1mode	Peaks 🔻	F1 output data mode		
	Imode	Integral 👻	Fitted intensity meaning		
	Scale	Linear -	Scaling		
	LWF	1	Line width factor	Ŧ	

#### Follow the indications shown in the red boxes!

1 Gyakorlat_Pr_190923 3 1 C:\Users\FANNI\Documents\NMR\spektrum						
Spectrum Pr	ocPars AcquPars Ti	tle PulseProg Peaks	Integrals Sample Structure Plot Fid			
PG <sup>L</sup>	) I, 🏓 🛄 🛡 🔗	}				
General First	General			·		
Second	Method	exponential -	Processing method	=		
Third	ExpVar	Gradient -	Variable parameter			
Baseline	Xlist	difflist	Variable parameter values file name			
Contain	Nstart	0	Start of input points			
	Ndata	8	Number of input points (TD)			
	Maxiter	100	Maximum number of iterations			
	EPS	1	Tolerance			
	Nexp	1	Number of components to fit			
	Noise	24351.00	Noise level (S_DEV)			
	PC	4	Noise sensitivity factor			
	SpiSup	1	Spike suppression factor			
	F1mode	Peaks -	F1 output data mode			
	Imode	Integral -	Fitted intensity meaning			
	Scale	Linear 🔹	Scaling			
	LWF	1	Line width factor	-		

Change scale instead of Linear to Logarithmic!

#### The following entries should be updated

![](_page_37_Picture_1.jpeg)

GAMMA: 5173 Hz/G

If all is set, then first click on button 1 - performs the calculation, then go to Spectrum (2) and look at the result Zoom on the SDS signals and read the log*D* value with positioning the cursor in the middle of the signals. On the upper left corner you can read the precise value of the y axis position ("Row")

	3 Gyakorlat_Pr_1	190923 3 1 C:\Users\FANN	II\Documents\NMR\spektrur	m 💼	• 💌
-	Spectrum Pr	ocPars AcquPars Ti	tle PulseProg Peaks	Integrals Sample Structure Plot Fid	
/	P G 🗓	<u>I, I, H</u> []] V 🔗			
1///////-	General	SpiSup	1	Spike suppression factor	*
' /	First	F1mode	Peaks -	F1 output data mode	
	Second	Imode	Integral 🔹	Fitted intensity meaning	
	Third	Scale	Logarithmic -	Scaling	
	Baseline	LWF	1	Line width factor	
	Conun	DISPmin	-10	Lower display limit	=
2		DISPmax	-8	Upper display limit	
		Npars	7	Number of parameters	
		Nvar	2	Number of parameters to fit	
		Gamma [Hz/G]	5173.00000	Gamma	
		Grad [G/cm]	0	Diffusion gradient	
		Gdist [ms]	150.00000	Gradient distance, big delta	
		Glen [ms]	2.00000	Gradient length, little delta	
	First component				
		l1vary	Yes 🔹	Fit intensity?	
		11	100000000	Intensity	
		14	A	and the second	Ŧ

An example can be found on slide 35.

#### The report should contain:

- Your names and the date
- Title, and goal of the practical
- Short theoretical description (SDS, micelles, cmc and ways of determination, DOSY, Stejskal-Tanner-equation)
- Steps of the measurement
- Stock solutions, sample preparations
- diffusion parameters  $(\delta, \Delta, \gamma, p1, o1)$
- eddosy 2D lg D
- 1D spectra for each solution, discussion
- Show a Stejskal-Tanner fitting from your data
- Results in table (c, chemical shift (ppm), D, average, standard deviation), 1D (perform it for both area and intensity), results from 2D evaluation
- Concentration-*D* curves, calculated cmc
- Check the literature values, comment the discrepancies
- Comments

Send the report in 1 week, email: abodor@caesar.elte.hu

### **Recommended literature**

https://www.cis.rit.edu/htbooks/nmr/inside.htm

#### The Basics of NMR

Joseph P. Hornak, Ph.D.

![](_page_40_Picture_4.jpeg)

Copyright © 1997-2019 J.P. Hornak. All Rights Reserved.